

## Product Datasheet

# Human High Sensitivity GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit (orb1945903)

<b>Description</b>	Human High Sensitivity GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.78-50pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma and Other biological fluids,Sample Volume= 100µL
<b>Sensitivity</b>	0.47 pg/mL
<b>Expiration Date</b>	Please enquire.