

## Product Datasheet

# Human High Sensitivity PGF(Placental Growth Factor) ELISA Kit (orb1945914)

<b>Description</b>	Human High Sensitivity PGF(Placental Growth Factor) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56-100pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Tissue homogenate and Other biological samples□Sample Volume=100μL
<b>Sensitivity</b>	0.94 pg/mL
<b>Expiration Date</b>	Please enquire.