



## Product Datasheet

# SARS-CoV-2 Wild Type Variant Recombinant Spike Trimer His Tag (orb1945933)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | SARS-CoV-2 Wild Type Variant Recombinant Spike Trimer His Tag   |
| <b>Reactivity</b>          | Other   |
| <b>Tested Applications</b> | ELISA   |
| <b>Preservatives</b>       | 1X PBS, 0.025% NaN3   |
| <b>Concentration</b>       | 0.5 mg/ml   |
| <b>Storage</b>             | 2-8°C   |
| <b>Tag</b>                 | Unconjugated  |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | SARS-CoV-2  |
| <b>Dilution Range</b>      | Control antigen for SARS-CoV-2 (COVID-19) diagnostic assays. Contains C-terminal His tag used for purification. |
| <b>Expiration Date</b>     | 6 months from date of receipt.  |



Concentration-response curves for bind...



SARS-CoV-2 Spike Trimer, His Tag (Wild T...