


Product Datasheet


SARS-CoV-2 BA.2.75 Variant Recombinant Spike RBD His Tag-20 (orb1945947)

Description	SARS-CoV-2 BA.2.75 Variant Recombinant Spike RBD His Tag-20
Reactivity	Other
Tested Applications	ELISA
Preservatives	1X PBS, 0.025% NaN3
Concentration	0.5 mg/ml
Storage	2-8°C
Tag	His tag
Note	For research use only
Isotype	SARS-CoV-2
Dilution Range	Control antigen for SARS-CoV-2 (COVID-19) diagnostic assays. Contains C-terminal His tag used for purification.
Expiration Date	6 months from date of receipt.



Concentration-response curves for bind...

SARS-CoV-2 BA.2.75 Spike Trimer Protein



SARS-CoV-2 Spike Trimer, His Tag (BA.2.7...