

## Product Datasheet

**SARS-CoV2 LxC1-G4 human mAb Antibody  
(orb1945969)**

**Description**

SARS-CoV2 LxC1-G4 human mAb Antibody

**Reactivity**

Other

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, FA

**Immunogen**

A recombinant SARS-CoV-2 spike protein containing RBD.

**Preservatives**

1X PBS

**Concentration**

0.5 mg/ml

**Storage**

-20°C

**Note**

For research use only

**Isotype**

Rabbit IgGk

**Clonality**


Monoclonal

**Dilution Range**


Experimental conditions may be optimized.

**Expiration Date**


12 months from date of receipt.



Spike protein trimers were coated at 2 µg/ml



Spike protein trimers were coated at 2 µg/ml



Spike protein trimers were coated at 2 µg/ml