

Product Datasheet

SARS-CoV2 LxC1-E7 human mAb Antibody (orb1945970)


Description

SARS-CoV2 LxC1-E7 human mAb Antibody

Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, FA
Immunogen	A recombinant SARS-CoV-2 spike protein containing RBD.
Preservatives	1X PBS
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Dilution Range	Experimental conditions may be optimized.
Expiration Date	12 months from date of receipt.



Spike protein trimers were coated at 2 µg/ml



Spike protein trimers were coated at 2 µg/ml



SARS-CoV2 LxC1-E7 human neutralizing mAb