
Product Datasheet

SARS-CoV-2 NP (84C4a) rabbit mAb Antibody (orb1945990)

Description

SARS-CoV-2 NP (84C4a) rabbit mAb Antibody

Reactivity

Other

Conjugation

Unconjugated

Tested Applications

ELISA

Immunogen

SARS-CoV-2 nucleoprotein (NP) specific peptide

Preservatives

1X PBS, 0.02% NaN3

Concentration

0.5 mg/ml

Storage

2-8°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality


Monoclonal

Uniprot ID
PODTC9
Dilution Range


Used for specific, high-sensitivity detection of SARS-CoV-2 nucleocapsid protein (NP) in immunoassay. Can be paired with other NP-specific clones for detection in sandwich ELISA format.

Expiration Date


12 months from date of receipt.



Western blot analysis of SARS-CoV2 and S...



Microtiter wells were coated with SARS-C...



A sandwich ELISA was performed using SAR...