

---


## Product Datasheet

### **SARS-CoV-2 NP (75G5a) rabbit mAb Biotin Conjugate Antibody (orb1945991)**


**Description**

SARS-CoV-2 NP (75G5a) rabbit mAb Biotin


<b>Reactivity</b>	Other
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	SARS-CoV-2 nucleoprotein (NP) specific peptide
<b>Preservatives</b>	1X PBS, 0.02% NaN3
<b>Concentration</b>	0.5 mg/ml
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgGk
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>PODTC9</b>
<b>Dilution Range</b>	Used for specific, high-sensitivity detection of SARS-CoV-2 nucleocapsid protein (NP) in immunoassay. Can be paired with other NP-specific clones for detection in sandwich ELISA format.
<b>Expiration Date</b>	12 months from date of receipt.



Microtiter wells were coated with SARS-CoV-2 NP...



Western blot analysis of SARS-CoV-2 NP and SARS-CoV-2...



A sandwich ELISA was performed using SARS-CoV-2 NP...