

Product Datasheet

Phospho-Jak3 (Tyr980/981) (E10) rabbit mAb FITC Conjugate Antibody (orb1945995)

Description Phospho-Jak3 (Tyr980/981) (E10) rabbit mAb FITC Conjugate Antibody

Reactivity Human

Conjugation FITC

Tested Applications FC

Immunogen A synthetic phospho-peptide corresponding to residues surrounding Tyr980/981 of human Jak3

Preservatives 1X PBS, 0.09% NaN₃, 0.2% rAlbumin

Storage 2-8°C

Note For research use only


Isotype Rabbit IgGk

Clonality Monoclonal

Uniprot ID [P52333](#)

Dilution Range For flow cytometric staining, the suggested use of this reagent is 5 µL per million cells or 5 µL per 100 µL of staining volume. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date 12 months from date of receipt.



Flow cytometric analysis of Ramos cells ...