



Product Datasheet

Phospho-NFKB p65 (Ser468) (B9) rabbit mAb Antibody (orb1945999)

Description	Phospho-NFKB p65 (Ser468) (B9) rabbit mAb Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Ser536 of human phospho-NFKB p65
Preservatives	1X PBS, 0.09% NaN3, 0.2% rAlbumin
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	Q04206
Dilution Range	For flow cytometric staining, the suggested use of this reagent is 5 µL per million cells or 5 µL per 100 µL of staining volume. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of A431 cells U...



Peptide blocking flow cytometric analysis...