


Product Datasheet

Phospho-PKC (pan) (gamma Thr514) (PF4) rabbit mAb FITC conjugate Antibody (orb1946006)

Description	Phospho-PKC (pan) (gamma Thr514) (PF4) rabbit mAb FITC conjugate Antibody
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Thr514 of human phospho PKC γ
Preservatives	1X PBS, 0.09% NaN ₃ , 0.2% rAlbumin
Storage	2-8°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	P05771, Q05655, P17252, P24723, Q04759, P05129, Q02156
Dilution Range	For flow cytometric staining, the suggested use of this reagent is 5 μ L per million cells or 5 μ L per 100 μ L of staining volume. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of HT1080 cells...