


## Product Datasheet

# Phospho-PKC (pan) (gamma Thr514) (PF4) rabbit mAb PE conjugate Antibody (orb1946007)

<b>Description</b>	Phospho-PKC (pan) (gamma Thr514) (PF4) rabbit mAb PE conjugate Antibody
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	PE
<b>Tested Applications</b>	FC
<b>Immunogen</b>	A synthetic phospho-peptide corresponding to residues surrounding Thr514 of human phospho PKC $\gamma$
<b>Preservatives</b>	1X PBS, 0.09% NaN <sub>3</sub> , 0.2% rAlbumin
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgGk
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P05771, Q05655, P17252, P24723, Q04759, P05129, Q02156</b>
<b>Dilution Range</b>	For flow cytometric staining, the suggested use of this reagent is 5 $\mu$ L per million cells or 5 $\mu$ L per 100 $\mu$ L of staining volume. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
<b>Expiration Date</b>	12 months from date of receipt.



Flow cytometric analysis of HT1080 cells...