


Product Datasheet


Phospho-SAPK/JNK (Thr183/Tyr185) (A11) rabbit mAb Antibody (orb1946073)

Description	Phospho-SAPK/JNK (Thr183/Tyr185) (A11) rabbit mAb Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Thr183/Tyr185 of human phospho SAPK/JNK
Preservatives	1X PBS, 0.02% NaN3, 50% Glycerol, 0.1% rAlbumin
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	P45983
Dilution Range	1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



APC	Fluor	Peptide	Block	Block	Block
APC	APC	None	None	None	None
APC	APC	None	None	None	None
APC	APC	None	None	None	None
APC	APC	None	None	None	None
APC	APC	None	None	None	None
APC	APC	None	None	None	None
APC	APC	None	None	None	None
APC	APC	None	None	None	None
APC	APC	None	None	None	None
APC	APC	None	None	None	None

Peptide blocking flow cytometric analysis...



Flow cytometric analysis of 293T cells S...