


Product Datasheet

Phospho-Shp1 (Tyr536) (2A7) rabbit mAb Antibody (orb1946077)

Description	Phospho-Shp1 (Tyr536) (2A7) rabbit mAb Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Tyr536 of human phospho Shp1
Preservatives	1X PBS, 0.02% NaN ₃ , 50% Glycerol, 0.1% rAlbumin
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	P29350
Dilution Range	1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of Ramos cells ...