


Product Datasheet


Phospho-PLK1 (Thr210) (C2) rabbit mAb Antibody (orb1946081)

Description	Phospho-PLK1 (Thr210) (C2) rabbit mAb Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Thr210 of human phospho-PLK1
Preservatives	1X PBS, 0.02% NaN3, 50% Glycerol, 0.1% rAlbumin
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	P53350
Dilution Range	1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



Fluorochrome	Phospho-PLK1 (T210)	Isotype
APC	100%	100%
APC	100%	100%
APC	100%	100%
APC	100%	100%
APC	100%	100%
APC	100%	100%
APC	100%	100%
APC	100%	100%
APC	100%	100%

Peptide blocking flow cytometric analysis...



Flow cytometric analysis of HeLa cells S...