
Product Datasheet

Phospho-Btk (Tyr551) (G12) rabbit mAb Antibody (orb1946087)

Description

Phospho-Btk (Tyr551) (G12) rabbit mAb Antibody

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA, FC

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Tyr551 of human phospho Btk

Preservatives

 1X PBS, 0.02% Na₃N, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality


Monoclonal

Uniprot ID
[Q06187](#)
Dilution Range


1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date


12 months from date of receipt.



Flow cytometric analysis of Ramos cells ...



Flow cytometric analysis of Ramos cells ...



Peptide ELISA using BtkY551-G12 (Cat. #o...