
Product Datasheet

**Phospho-Stat5 (Tyr694) (G11) rabbit mAb
Antibody (orb1946100)**

Description

Phospho-Stat5 (Tyr694) (G11) rabbit mAb

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

FC, WB

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Tyr694 of human phospho Stat5

Preservatives

 1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality


Monoclonal

Uniprot ID
P51692, P42229
Dilution Range


1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date

12 months from date of receipt.




Western blot analysis of K562 cell extra...



Fluor	Isotype	Peptide	Block	Median
FL1	IgG	None	0	100
FL1	IgG	None	100	100
FL1	IgG	Phospho	0	100
FL1	IgG	Phospho	100	100
FL1	IgG	None	0	100
FL1	IgG	None	100	100
FL1	IgG	None	0	100
FL1	IgG	None	100	100

Peptide blocking flow cytometric analysis...



Flow cytometric analysis of 3T3 cells se...