


Product Datasheet

Phospho-NFkB p65 (Ser536) (C8) rabbit mAb Antibody (orb1946111)

Description	Phospho-NFkB p65 (Ser536) (C8) rabbit mAb Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Ser536 of human phospho NFkB.
Preservatives	1X PBS, 0.02% NaN3, 50% Glycerol, 0.1% rAlbumin
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	Q04206
Dilution Range	1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell extra...