

---

## Product Datasheet

### Phospho-BLNK (Tyr84) (H4) rabbit mAb Antibody (orb1946129)

**Description**

Phospho-BLNK (Tyr84) (H4) rabbit mAb Antibody

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

FC

**Immunogen**

A synthetic phospho-peptide corresponding to residues surrounding Tyr84 of human phospho BLNK

**Preservatives**

 1X PBS, 0.02% Na<sub>3</sub>N, 50% Glycerol, 0.1% rAlbumin

**Concentration**

0.5 mg/ml

**Storage**

-20°C

**Note**

For research use only

**Isotype**

Rabbit IgGk

**Clonality**

Monoclonal

**Uniprot ID**
**Q8WV28**
**Dilution Range**

1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.


**Expiration Date**

12 months from date of receipt.




ISG	Treatment	Peptide Block	Median	SD
ISG1	IFNα	Isotype	8653	128
ISG2	IFNα	Isotype	128	128
ISG3	IFNα	Phos	108	128
ISG4	IFNα	Phos	262	128
ISG5	IFNα	Phos	2628	128
ISG6	IFNα	Phos	4192	128
ISG7	IFNα	Phos	5756	128

Peptide blocking flow cytometric analysis...



Flow cytometric analysis of Daudi cells



Flow cytometric analysis of Daudi cells