
Product Datasheet

Phospho-Stat5 (Tyr694) (B5) rabbit mAb Antibody (orb1946160)

Description

Phospho-Stat5 (Tyr694) (B5) rabbit mAb

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

FC

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Tyr694 of human phospho Stat5

Preservatives

1X PBS, 0.02% NaN3, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality


Monoclonal

Uniprot ID
P51692, P42229
Dilution Range


1.00 - 0.001 µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date


12 months from date of receipt.



Flow cytometric analysis of Ramos cells ...



Peptide blocking flow cytometric analysis...



Flow cytometric analysis of NIH3T3 cells...