
Product Datasheet

**Phospho-PLC γ 1 (Tyr783) (C4) rabbit mAb Antibody
(orb1946203)**

Description

Phospho-PLCγ1 (Tyr783) (C4) rabbit mAb


Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

FC

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Tyr783 of human phospho PLCγ1.

Preservatives

 1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality


Monoclonal

Uniprot ID
P19174
Dilution Range

1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.


Expiration Date

12 months from date of receipt.




SWELLID	Treatment	Median	SI
CA-APP	PK	3180	
CA-NM	PK	27480	
CA-D-05 (pp-PLC)	PK	25000	
CA-D-05 (pp-PLC)	Imatinib	2784	
CA-D-05 (pp-PLC)	Imatinib	699	

Peptide blockage flow cytometric analysis...



PLCγ1Y783-C4 recognizes basal phosphoryl...



Flow cytometric analysis of HeLa cells S...