


## Product Datasheet

# Phospho-KSR1 (Ser392) (3A4) rabbit mAb APC conjugate-10 Tests Antibody (orb1946214)

<b>Description</b>	Phospho-KSR1 (Ser392) (3A4) rabbit mAb APC conjugate-10 Tests Antibody
<b>Reactivity</b>	Human
<b>Conjugation</b>	APC
<b>Tested Applications</b>	FC
<b>Immunogen</b>	A synthetic phosphor-peptide corresponding to residues surrounding Ser392 of human phospho KSR1
<b>Preservatives</b>	1X PBS, 0.09% NaN <sub>3</sub> , 0.2% rAlbumin
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgGk
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>Q8IVT5</b>
<b>Dilution Range</b>	For flow cytometric staining, the suggested use of this reagent is 5 µL per million cells or 5 µL per 100 µL of staining volume. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
<b>Expiration Date</b>	12 months from date of receipt.



Flow cytometric analysis of HeLa cells u...