
Product Datasheet

**Phospho-Chk1 (Ser317) (F10) rabbit mAb Antibody
(orb1946237)**

Description

Phospho-Chk1 (Ser317) (F10) rabbit mAb


Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

FC

Applications
Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Ser317 of human phospho Chk1

Peptide blocking flow cytometric analysis...

Preservatives

1X PBS, 0.02% NaN3, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality


Monoclonal

Uniprot ID
O14757
Dilution Range


1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date

12 months from date of receipt.



Flow cytometric analysis of NIH3T3 cells...



Flow cytometric analysis of HeLa cells S...