


Product Datasheet

Isotype Control (G9) rabbit mAb FITC Conjugate Antibody (orb1946239)

Description	Isotype Control (G9) rabbit mAb FITC Conjugate Antibody
Reactivity	Other
Conjugation	FITC
Tested Applications	FC
Immunogen	Keyhole limpet hemocyanin carrier protein.
Preservatives	1X PBS, 0.09% NaN ₃ , 0.2% rAlbumin
Storage	2-8°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Dilution Range	Isotype control antibodies should be used at the same concentration as the test antibody.
Expiration Date	12 months from date of receipt.



Isotype Control FITC

Flow cytometric analysis of U937 cells u...