

Product Datasheet

Isotype Control (G9) rabbit mAb PE Conjugate (orb1946240)


Description	Isotype Control (G9) rabbit mAb PE Conjugate
Reactivity	Other
Conjugation	PE
Tested Applications	FC
Immunogen	Keyhole limpet hemocyanin carrier protein.
Preservatives	1X PBS, 0.09% NaN ₃ , 0.2% rAlbumin
Storage	2-8°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Clone Number	KLH-G9
Dilution Range	Isotype control antibodies should be used at the same concentration as the test antibody.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Flow cytometric analysis of U937 cells using rabbit (G9) mAb Isotype Control (PE Conjugate) for untreated (blue) or treated with IFNa (Red) compared to Phospho-Src (Tyr416) (C4) rabbit mAb (PE conjugate #orb1946292) untreated (black) or IFNa treated (green). Cat. #orb1946240.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)