
Product Datasheet

**Phospho-NDRG1 (Thr346) (F5) rabbit mAb
Antibody (orb1946264)**

Description

Phospho-NDRG1 (Thr346) (F5) rabbit mAb


Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

FC, WB

Applications
Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Thr346 of human phospho NDRG1

Preservatives

1X PBS, 0.02% NaN3, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality

Monoclonal

Uniprot ID
[Q92597](#)
Dilution Range

1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.


Expiration Date

12 months from date of receipt.




IS	IS	IS	IS	IS	IS	IS	IS	IS	IS
101	101	101	101	101	101	101	101	101	101
102	102	102	102	102	102	102	102	102	102
103	103	103	103	103	103	103	103	103	103
104	104	104	104	104	104	104	104	104	104
105	105	105	105	105	105	105	105	105	105
106	106	106	106	106	106	106	106	106	106
107	107	107	107	107	107	107	107	107	107
108	108	108	108	108	108	108	108	108	108
109	109	109	109	109	109	109	109	109	109
110	110	110	110	110	110	110	110	110	110
111	111	111	111	111	111	111	111	111	111
112	112	112	112	112	112	112	112	112	112
113	113	113	113	113	113	113	113	113	113
114	114	114	114	114	114	114	114	114	114
115	115	115	115	115	115	115	115	115	115
116	116	116	116	116	116	116	116	116	116
117	117	117	117	117	117	117	117	117	117
118	118	118	118	118	118	118	118	118	118
119	119	119	119	119	119	119	119	119	119
120	120	120	120	120	120	120	120	120	120

Peptide blocking flow cytometric analysis...



Flow cytometric analysis of THP1 cells S...



Western blot analysis of C2C12 cell extr...