
Product Datasheet

**Phospho-EGFR (Tyr1068) (E5) rabbit mAb Antibody
(orb1946283)**

Description

Phospho-EGFR (Tyr1068) (E5) rabbit mAb

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Tyr1068 of human phospho EGFR.

Preservatives

 1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality


Monoclonal

Uniprot ID
P00533
Dilution Range


1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date


12 months from date of receipt.



Flow cytometric analysis of C6 cells sec...



Flow cytometric analysis of HeLa cells S...



Flow cytometric analysis of 3T3 cells se...