

Product Datasheet

Phospho-Zap70 (Tyr319)/Syk (Tyr352) (A3) rabbit mAb SureLight 488 conjugate Antibody (orb1946284)

Description Phospho-Zap70 (Tyr319)/Syk (Tyr352) (A3) rabbit mAb SureLight 488 conjugate Antibody

Reactivity Human, Mouse

Conjugation SureLight 488

Tested Applications FC

Immunogen A synthetic phospho-peptide corresponding to residues surrounding Tyr319/Tyr352 of human phospho Zap70/Syk.

Preservatives 1X PBS, 0.09% NaN₃, 0.2% rAlbumin

Storage 2-8°C

Note For research use only


Isotype Rabbit IgGk

Clonality Monoclonal

Uniprot ID [P43405](#), [P43403](#)

Dilution Range For flow cytometric staining, the suggested use of this reagent is 5 µL per million cells or 5 µL per 100 µL of staining volume. It is recommended that the reagent be titrated for optimal performance for each application.

Expiration Date 12 months from date of receipt.



Flow cytometric analysis of Jurkat cells...