

Product Datasheet

Phospho-Zap70 (Tyr319)/Syk (Tyr352) (A3) rabbit mAb Antibody (orb1946288)

Description Phospho-Zap70 (Tyr319)/Syk (Tyr352) (A3) rabbit mAb Antibody

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications FC

Immunogen A synthetic phospho-peptide corresponding to residues surrounding Tyr319/Tyr352 of human phospho Zap70/Syk.

Preservatives 1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration 0.5 mg/ml

Storage -20°C

Note For research use only


Isotype Rabbit IgGk

Clonality Monoclonal


Uniprot ID **P43405, P43403**

Dilution Range 1µg/mL – 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.


Expiration Date 12 months from date of receipt.



ZapY319-A3 recognizes basal phosphorylat...



Flow cytometric analysis of Jurkat cells...



Flow cytometric analysis of Jurkat cells...