


Product Datasheet


**Phospho-PI3 Kinase p85 (Tyr458)/p55 (Tyr199)
(1A11) rabbit mAb Antibody (orb1946298)**

Description


Phospho-PI3 Kinase p85 (Tyr458)/p55



Western
blot
analysis of
NIH3T3
cell ext...



Flow
cytometric
analysis of
NIH3T3
cells...



Peptide
blocking
flow
cytometric
analysi...

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Tyr458 of human phospho PI3K p85

Preservatives

1X PBS, 0.02% NaN3, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality

Monoclonal

Uniprot ID

P27986, **O00459**, **Q92569**

Dilution Range

1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date

12 months from date of receipt.