

Product Datasheet

Phospho-p53 (Ser15) (1C11) rabbit mAb Antibody (orb1946331)

Description

Phospho-p53 (Ser15) (1C11) rabbit mAb

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

FC, WB

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Ser15 of human phospho p53.

Preservatives

1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality

Monoclonal

Uniprot ID


P04637

Dilution Range


1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date


12 months from date of receipt.



Western blot analysis of L929 cell extra...



Flow cytometric analysis of C2C12 cells ...



Flow cytometric analysis of SK-N-MC cell...