

Product Datasheet

**Phospho-PLC γ 2 (Tyr759) (G3) rabbit mAb Antibody
(orb1946341)**

Description

 Phospho-PLC γ 2 (Tyr759) (G3) rabbit mAb

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

FC, WB

Immunogen

 A synthetic phospho-peptide corresponding to residues surrounding Tyr759 of human phospho PLC γ 2.

Preservatives

 1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality

Monoclonal


Uniprot ID
P16885
Dilution Range

 1 μ g/mL - 0.001 μ g/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.


Expiration Date

12 months from date of receipt.

Western blot analysis of Ramos cell extr...



Flow cytometric analysis of NIH3T3 cells...



Flow cytometric analysis of Ramos cells ...