


Product Datasheet

Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204) rabbit mAb SureLight 488 conjugate Antibody (orb1946366)

Description	Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204) rabbit mAb SureLight 488 conjugate Antibody
Reactivity	Human, Mouse
Conjugation	SureLight 488
Tested Applications	FC
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Thr202/Tyr204 of human phospho Erk1/2.
Preservatives	1X PBS, 0.09% NaN3, 0.2% rAlbumin
Storage	2-8°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	P28482 , P27361
Dilution Range	For flow cytometric staining, the suggested use of this reagent is 5 µL per million cells or 5 µL per 100 µL of staining volume. It is recommended that the reagent be titrated for optimal performance for each application.
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of Jurkat cells...