

## Product Datasheet

### Chicken GHRL(Ghrelin) ELISA Kit (orb1946456)

<b>Description</b>	Chicken GHRL(Ghrelin) ELISA Kit
<b>Reactivity</b>	Gallus
<b>Range</b>	156.25-10000 pg/mL
<b>Target</b>	GHRL
<b>Note</b>	For research use only
<b>Application notes</b>	<p>This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Ghrelin(GHRL) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Ghrelin(GHRL). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm <math>\pm</math> 10nm. The concentration of Ghrelin(GHRL) in the samples is then determined by comparing the OD of the samples to the standard curve.</p>
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Sensitivity</b>	46.4 pg/mL
<b>Expiration Date</b>	Please enquire.

