

## Product Datasheet

### Rabbit PRL(Prolactin) ELISA Kit (orb1946458)

<b>Description</b>	Rabbit PRL(Prolactin) ELISA Kit
<b>Reactivity</b>	Rabbit
<b>Range</b>	3.13-200 ng/mL
<b>Target</b>	PRL
<b>Note</b>	For research use only
<b>Application notes</b>	<p>This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Prolactin(PRL) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Prolactin(PRL). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Prolactin(PRL) in the samples is then determined by comparing the OD of the samples to the standard curve.</p>
<b>Sample Types</b>	Serum, plasma and other biological fluids
<b>Sensitivity</b>	0.92 ng/mL
<b>Expiration Date</b>	Please enquire.

