

## Product Datasheet

### Chicken GC(Glucagon) ELISA Kit (orb1946475)

<b>Description</b>	Chicken GC(Glucagon) ELISA Kit
<b>Reactivity</b>	Gallus
<b>Range</b>	23.44-1500 pg/mL
<b>Target</b>	GC
<b>Note</b>	For research use only
<b>Application notes</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Glucagon(GC) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Glucagon(GC). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Glucagon(GC) in the samples is then determined by comparing the OD of the samples to the standard curve.
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Sensitivity</b>	7.59 pg/mL
<b>Expiration Date</b>	Please enquire.

