

Product Datasheet

Guinea Porcine GT(Gastrin) ELISA Kit [Out of stock] (orb1946491)

Description	Guinea Porcine GT(Gastrin) ELISA Kit [Out of stock]
Reactivity	Guinea pig
Range	15.63-1000 pg/mL
Target	GT
Note	For research use only
Application notes	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Gastrin(GT) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Gastrin(GT). Next,Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Gastrin(GT) in the samples is then determined by comparing the OD of the samples to the standard curve.
Sample Types	Serum, plasma and other biological fluids.
Sensitivity	4.35 pg/mL
Expiration Date	Please enquire.

