

Product Datasheet

Sheep Cor(Cortisol) ELISA Kit (orb1946593)

Description	Sheep Cor(Cortisol) ELISA Kit
Reactivity	Sheep
Range	3.13-200 ng/mL
Target	Cor
Note	For research use only
Application notes	<p>This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Cortisol(Cor) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cortisol(Cor). Next,Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Cortisol(Cor) in the samples is then determined by comparing the OD of the samples to the standard curve.</p>
Sample Types	serum, plasma, urine and other biological fluids
Sensitivity	0.93 ng/mL
Expiration Date	Please enquire.

