

## Product Datasheet

### Rabbit Cr(Calretinin) ELISA Kit (orb1946619)

<b>Description</b>	Rabbit Cr(Calretinin) ELISA Kit
<b>Reactivity</b>	Rabbit
<b>Range</b>	1.25-80 nmol/ml
<b>Target</b>	CR
<b>Note</b>	For research use only
<b>Application notes</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Calretinin(Cr) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Calretinin(Cr). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Calretinin(Cr) in the samples is then determined by comparing the OD of the samples to the standard curve.
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Sensitivity</b>	0.75 nmol/ml
<b>Expiration Date</b>	Please enquire.

