

Product Datasheet

Cattle TRH(Thyrotropin Releasing Hormone) ELISA Kit (orb1946636)

Description	Cattle TRH(Thyrotropin Releasing Hormone) ELISA Kit
Reactivity	Bovine
Range	7.82-500 pg/mL
Target	TRH
Note	For research use only
Application notes	<p>This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Thyrotropin Releasing Hormone(Cattle TRH) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Thyrotropin Releasing Hormone(Cattle TRH). Next,Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Thyrotropin Releasing Hormone(Cattle TRH) in the samples is then determined by comparing the OD of the samples to the standard curve.</p>
Sample Types	serum, plasma and other biological fluids
Sensitivity	2.63 pg/mL
Expiration Date	Please enquire.

