

## Product Datasheet

### DC(Deoxycholate) ELISA Kit (orb1946735)

<b>Description</b>	DC(Deoxycholate) ELISA Kit
<b>Range</b>	15.63-1000 ng/mL
<b>Target</b>	DC
<b>Note</b>	For research use only
<b>Application notes</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Deoxycholate(DC) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Deoxycholate(DC). Next,Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of Deoxycholate(DC) in the samples is then determined by comparing the OD of the samples to the standard curve.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Sensitivity</b>	4.99 ng/mL
<b>Expiration Date</b>	Please enquire.

