


---

## Product Datasheet

### Mouse TRH(Thyrotropin Releasing Hormone) ELISA Kit (orb1946772)

**Description**

Mouse TRH(Thyrotropin Releasing


**Reactivity**

Mouse

**Range**

7.82-500pg/mL

**Target**

TRH

**Note**

For research use only

**Application notes**

This ELISA kit uses the Competitive-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with Mouse TRH. During the reaction, Mouse TRH in the sample or standard competes with a fixed amount of Mouse TRH on the solid phase supporter for sites on the Biotinylated Detection Ab specific to Mouse TRH. Excess conjugate and unbound sample or standard are washed away, and Avidin-Horseradish Peroxidase (HRP) conjugate are added to each micro plate well and incubated. Then a TMB substrate solution is added to each well. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns from blue to yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of  $450 \text{ nm} \pm 2 \text{ nm}$ . The concentration of Mouse TRH in tested samples can be calculated by comparing the OD of the samples to the standard curve.

**Sample Types**

serum, plasma, Tissue homogenate and Other biological samples

**Uniprot ID**
[Q62361](#)
**Sensitivity**

4.69 pg/mL

**Expiration Date**

Please enquire.