

Product Datasheet

Rat T4-Ab(Thyroxine Antibody) ELISA Kit (orb1947413)

Description	Rat T4-Ab(Thyroxine Antibody) ELISA Kit
Reactivity	Rat
Range	1.57-100pmol/mL
Target	T4-Ab
Note	For research use only
Application notes	This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Rat T4-Ab. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Rat T4-Ab and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Rat T4-Ab, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm $\pm$ 2 nm. The OD value is proportional to the concentration of Rat T4-Ab. You can calculate the concentration of Rat T4-Ab in the samples by comparing the OD of the samples to the standard curve.
Sample Types	serum, plasma, Tissue homogenate and Other biological samples
Sensitivity	0.94 pmol/mL
Expiration Date	Please enquire.

