



Product Datasheet

SCN9A Antibody (orb194958)

Description	SCN9A Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q15858
Expiration Date	12 months from date of receipt.



Western blot analysis of RAW cells using...



Immunohistochemistry analysis of human I...