

Product Datasheet

EasyStep Mouse T4(Thyroxine) ELISA Kit (orb1950005)

Description	EasyStep Mouse T4(Thyroxine) ELISA Kit
Reactivity	Mouse
Range	9.38-600 ng/mL
Note	For research use only
Application notes	<p>This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Mouse T4 antigen, and the Mouse T4 standard plate wells that pre-coated using protein-related techniques are provided separately. Standard/Sample Diluent Buffer or samples are added to the appropriate microtiter plate wells ,then added a HRP-conjugated antibody specific to Mouse T4 to each microplate well and incubated . After TMB substrate solution is added, the enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm \pm10nm. The concentration of Mouse T4 in the samples is then determined by comparing the OD of the samples to the standard curve.</p>
Sample Types	Serum, plasma
Sensitivity	7 ng/mL
Expiration Date	Please enquire.