

Product Datasheet

CSDE1 Rabbit mAb Antibody (orb1950251)

Description	RNA-binding protein. Required for internal initiation of translation of human rhinovirus RNA. May be...
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	A synthesized peptide derived from human CSDE1
Form/Appearance	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	O75534
Dilution Range	WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50
Expiration Date	12 months from date of receipt.