

## Product Datasheet

# Retinoid X Receptor alpha Rabbit mAb Antibody (orb1950635)

<b>Description</b>	Receptor for retinoic acid. Retinoic acid receptors bind as heterodimers to their target response el...
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IP, WB
<b>Immunogen</b>	A synthesized peptide derived from human Retinoid X Receptor alpha
<b>Form/Appearance</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P19793</b>
<b>Dilution Range</b>	WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:50
<b>Expiration Date</b>	12 months from date of receipt.