

Product Datasheet

KMT1B / SUV39H2 Rabbit mAb Antibody (orb1950734)

| | |
|----------------------------|--|
| Description | KMT1B / SUV39H2 Rabbit mAb Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC, IP, WB |
| Storage | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | Q9H511 |
| Expiration Date | 12 months from date of receipt. |