

Product Datasheet

Collagen, type VII antibody (orb19514)

Description	Goat polyclonal to Collagen, type VII
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Collagen, type VII
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This product does not work on cryostat sections of Normal Human Skin.ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 250kDa band observed in rat testis lysates (calculated MW of 159kDa according to Rat NP_001100328.1). Recommended concentration: 0.1-1.5 µg/ml.Experiment Notes: A customer found this antibody unsuccessful on cryosections of Human Skin.
Protein Sequence	TVQYSDDPRTEF
Clonality	Polyclonal
MW	295
Entrez	1294
Dilution Range	ELISA: 1:32000, WB: 0.1-1.5 µg/ml
Expiration Date	12 months from date of receipt.



orb19514
(0.5ug/ml)
staining of rat test...