

Product Datasheet

GFP Rabbit mAb(HRP conjugated) Antibody (orb1951994)

Description	GFP Rabbit mAb(HRP conjugated) Antibody
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q963I9
Expiration Date	12 months from date of receipt.